

# TOP GLOVE

TOP QUALITY, TOP EFFICIENCY

THE WORLD'S LARGEST MANUFACTURER OF GLOVES

## DOUBLE DONNING SURGICAL GLOVE



### Faster Identification

Green colour inner glove facilitate earlier detection when the outer glove start to fail



### Functional Benefits:

- Provide double reduction on risk of patients' blood exposure
- Provide double assurance to avoid cross-infection between users and patients
- Provide double awareness to user on any perforation of glove by colour identification

# DOUBLE DONNING SURGICAL GLOVE

## Quality Standards

- Conforms to ASTM D3577 Standards
- Manufactured under QSR (GMP), ISO9001:2015 Quality Management System and ISO 13485:2016
- Sterilized by Gamma Ray Irradiation
- Bioburden and sterility tested
- Low Protein reducing potential allergic reactions



## Glove Sizes

- 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0
- Marked in the check box on the shipping carton

## Product Specifications

- Type** : Latex surgical glove, powder free, sterile
- Material** : Natural high grade rubber latex
- Colour** : Natural and Green
- Design & Features** : Hand specific, curved fingers, palm textured, beaded cuff
- Storage** : The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight
- Shelf-life** : The gloves shall have shelf-life of 5 years from the date of manufacture



## Physical Dimensions

Dimensions	STANDARDS	
	Top Glove	ASTM D3577
<b>Length (mm)</b>	Min 280	Min 245 ( size 5.5) Min 265 ( size 6-9)
<b>Width (mm)</b>		
5.5	71 ± 5	70 ± 6
6.0	77 ± 5	76 ± 6
6.5	83 ± 5	83 ± 6
7.0	89 ± 5	89 ± 6
7.5	95 ± 5	95 ± 6
8.0	102 ± 6	102 ± 6
8.5	108 ± 6	108 ± 6
9.0	115 ± 5	114 ± 6
<b>Thickness - Single Wall (mm)</b>		
• Fingers	Min 0.10	Min 0.10
• Palm	Min 0.10	Min 0.10
• Cuff	Min 0.10	Min 0.10

## Physical Properties

Property	ASTM D3577
<b>Elongation at break (%)</b>	
• Before Aging	Min 750
• After Aging	Min 560
<b>Tensile Strength (MPa)</b>	
• Before Aging	Min 24
• After Aging	Min 18

If you have any enquiries on our products, please contact:

### TOP GLOVE SALES & CORPORATE OFFICE

Level 21, Top Glove Tower,  
16, Persiaran Setia Dagang,  
Setia Alam, Seksyen U13,  
40170 Shah Alam,  
Selangor D.E, Malaysia.

Tel : +603-3362 3098  
Fax : +603-3362 5096 (Sales)  
Email : sales@topglove.com.my



Scan for TG Website



<https://www.facebook.com/mytopglove>



[www.topglove.com](http://www.topglove.com)